



Be Right™



Syringe filter, non-sterile PDVF, 33 mm, 0.45 um, 50/pk

Nr. produs: 2513603
RON Preț (fără TVA): Contact
Expediere în 2 săptămâni

Hach syringe filters offer an efficient way to remove particulates from samples prior to analysis. Filters feature 20% more surface area than standard 25 mm filters, which increases flow rates and capacity. PVDF (polyvinylidene fluoride) filters are ideal for solutions containing proteins. Standard Female Luer-Lok inlet with standard male slip Luer outlet. Hydrophilic. Diameter: 33mm. Pore Size: 0.45um. 50/PK.

Specificații

Cantitate:	50 /pk
Diametru:	33 mm
Dimensiune por:	0.45 μm
Material:	PVDF
Unitate:	50/pk