



Test Kit, AccuVac® Color Disc, Ozone, 0.05-1.50 mg/L

Nr. produs:

2518050

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



This ozone test kit uses the fast, accurate visual method for determining ozone concentrations. Use to control disinfectant levels in water, wastewater, and industrial water. Kit contains 25 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

AccuVac® vials mean no mess, no measuring

Indigo trisulfonate method masks interferences from chlorine and other oxidizing agents and makes ozone analysis easier

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Specificații

Domeniu de măsurare:	0 - 1.50 mg/L O ₃
Interval de măsurare standard:	0 - 1.50 mg/L O ₃
Model:	NA
Model test kit:	AccuVac® Kit
Numar teste:	25
Nume metoda:	Color disc/ Indigo trisulfonate
Parameter:	Ozone
Parametru:	Ozone - As O ₃
Platforma:	Accuvac
Smallest increment:	0.05
Tipul carcasei:	G