



## Bottle rack, plastic-coated, 20 compartments

**Nr. produs:** 2551801

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak® sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

---

### Specificații

Cantitate:	1
Caracteristici speciale:	Rack, 20 Bottle Compartments