



Sample container, 120 mL, pk/250

Nr. produs:	2552212
RON Preț (fără TVA):	Contact
Expediere în 2 săptămâni	

Plastic, leak-proof. Corning Snap Seal offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with patented sealing mechanism so sample won't leak.

Specificații

Cantitate:	250 /pk
Capacitate:	120 mL - 250/pk
Diametru exterior :	45 mm (approx)