



Test kit, chloride, model CDS-DT

Nr. produs: 2580600
RON Preț (fără TVA): Contact
 Nu este afișată data expirării

Produse periculoase



The single-parameter Model CDS-DT uses the Hach Digital Titrator with Silver Nitrate Titrant Cartridges to determine chloride concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-10,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specificații

Domeniu de măsurare:	10 - 10,000 mg/L Cl ⁻
Domeniu de masurare :	to
Model test kit:	CDS-DT
Numar teste:	100
Nume metoda:	Digital Titrator/ Silver nitrate
Parametru:	Chloride - As Cl ⁻ , 10-10,000 mg/L Cl ⁻
Platforma:	Titrator digital
Smallest increment:	10
Tipul carcasei:	F

Accesorii necesare

- Pipetă, volum variabil, 1,0-5,0 mL (Item BBP065)