



## Standard solution, chromium, hexavalent, 5 mg/L as Cr<sup>6+</sup> (NIST), pk/20 - 2 mL PourRite® ampule

Nr. produs: 2605620

RON Preț (fără TVA): Contact

Nu este afișată data expirării

For accuracy checks in chromium determinations. Pack of 20 - 2 mL PourRite(TM) Ampules.

### Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

### Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

### Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

---

## Specificații

Cantitate:	2 mL
Concentrație:	5 mg/L
Descriere:	Chromium Hexavalent Standard Solution as Cr <sup>6+</sup> (HIST)
Platforma:	Solu#ie