



## Pitcher, graduated, 2000 cc

**Nr. produs:** 2612854  
**RON Preț (fără TVA):** Contact  
Nu este afișată data expedierii

Chemical-resistant, shatter-proof, tall-form pitchers with well-designed spouts, and large handles. Metric graduations are molded onto outside wall of each pitcher. Polypropylene, steam autoclavable at 250°F (121°C).

---

### Specificații

Adancime:	8.1 Inches
Cantitate:	1 each
Capacitate:	2000 mL
Graduations:	100 mL