



## Cylinder, graduated, 1000 mL, w/ handles

Nr. produs: 2612953

RON Preț (fără TVA): Contact

Nu este afișată data expirării

One-piece polypropylene construction with a molded in handle top ensure a secure, solid grip and superior control for safer, more convenient pouring. The “Total Contained, Total Delivered” double scale is accurate and easy to read (cylinder will not create a meniscus), while the hexagonal base prevents tipping or rolling. Hand opening is 89x 29 mm (3 x 11/8"). Steam autoclavable upto 121°C (250°F).

---

### Specificații

Cantitate:	1 each
Capacitate:	1000 mL
Certificări:	None
Material corp:	Polypropylene
Metoda de măsurare:	To Deliver/To Contain
Tip:	Dublu metric