



## Bottle, wash, diamond alert, 500 mL, Variety 6/pk

Nr. produs: 2613049  
RON Preț (fără TVA): Contact  
Expediere în 2 săptămâni

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

---

### Specificații

Cantitate:	6 pk
Capacitate:	500 mL
Culoare:	Transparent
Material:	LDPE
Stem:	Blue Cap/Stem
Unitate:	6 pk