



Filter, membrane, pore size 0.2 µm, diameter 47 mm, 200/pk (Gelman)

Nr. produs: 2618800

RON Preț (fără TVA): Contact

Nu este afișată data expirării

For isolation and enumeration of stressed organisms (*Pseudomonas* species). High retention and flow rate. Certified for use in membrane filtration analysis. White polyethersulfone media. STerile and individually packaged. Pore size: 0.2 µm. Diameter: 47mm.

Specificații

Cantitate:	200 /pk
Diametru:	47 mm
Dimensiune por:	0.2 µm
Filter Type:	Filter
Material:	Polyethersulfone