



Sample vial, 120 mL w/ thiosulfate & custody tie, 200/pk

Nr. produs: 2634001
RON Preț (fără TVA): Contact
Expediere în 2 săptămâni

Sample Vial with Thiosulfate and Custody Tie

Sample Vial with Thiosulfate and Custody Tie

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene.

Specificații

Cantitate:	200 /pk
Capacitate:	fill line at 100 mL (+/-2.5 mL)
Caracteristici speciale:	Sample Vial, with Thiosulfate
Diametru:	43mm
Lungime:	93mm
Unitate:	200 /pk