



Sample vial, 120 mL w/out thiosulfate, 200/pk

Nr. produs: 2634002
RON Preț (fără TVA): Contact
Nu este afișată data expedierii

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene.

Specificații

Cantitate:	200/pk
Capacitate:	Fill line at 100 mL ($\pm$ 2.5 mL)
Caracteristici speciale:	Sample Vial, without Thiosulfate
Diametru:	43 mm
Lungime:	93 mm