



## Flask, volumetric, class A, certified, 1000 mL

Nr. produs: 2636653

### ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.

Individually serialized glass flask supplied with Certificate of Graduation Accuracy and traceable to NIST Standards. Class A tolerances meet requirements of ASTM E438 and ASTM E288. Graduation ring on neck of flask. Supplied with ground glass stopper. (Stopper size 22).

---

### Specificații

Capacitate:	1000 mL
Certificări:	Yes
Material corp:	Class A Glass
Unitate:	1 each