



Test kit, RapidSilver™ visual

Nr. produs: 2674500

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.

Test Kit Silver, Model RapidSilver™. Range 0-50 ppb. Silver As Ag Visual Comparison.

Provides fast, accurate information on silver discharges

Used by surface finishers, photographic film manufacturers and processors, and others to prevent pollution and economic loss

Kit directions also include digestion and dilution procedures needed for testing silver strike bath samples

Specificații

Domeniu de măsurare:	5 - 50 ppb Ag
Interval de măsurare standard:	5 - 50 ppb Ag
Model:	NA
Model test kit:	RapidSilver™
Numar teste:	100
Nume metoda:	Visual Comparison
Parameter:	Silver
Parametru:	Silver-As Ag
Smallest increment:	5 µg/L
Tipul carcasei:	F