



## Test kit, hypochlorite, model CN-HR

**Nr. produs:**

**2687200**

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Hypochlorite Kit, Model CN-HR - Drop Count

Hypochlorite Kit, Model CN-HR - Drop Count

Utilizing the simple, economical drop count titration method, this kit is used to determine high range hypochlorite solution concentrations (bleach solutions). Range: 5-15%, contains 100 tests.

**Simple, economical titration method**

**Kit contains all required reagents, solutions, and glassware**

### Specificații

|                      |                                                      |
|----------------------|------------------------------------------------------|
| Domeniu de măsurare: | 50 - 150 g/L $\text{Cl}_2$ (5-15%)                   |
| Model test kit:      | CN-HR                                                |
| Numar teste:         | 100                                                  |
| Nume metoda:         | Drop count titration/<br>Thiosulfate                 |
| Parametru:           | Hypochlorite, high range - (Bleach) As $\text{Cl}_2$ |
| Platforma:           | Drop Count Titration                                 |
| Smallest increment:  | 0.5%                                                 |
| Tipul carcasei:      | D                                                    |