



SpecCheck Secondary Gel Standard Kit - Fluoride, 0-2.0 mg/L F

Nr. produs: 2712500

RON Preț (fără TVA): Contact
Expediere în 2 săptămâni

Solid, ready for use Gel-Standards for instrument verification in the lab and in the field.

Set contains four vials filled with gels that simulate the test colour at various concentrations.

Specificații

Concentrație:	0 – 2,00 mg/L F
Instrument:	All except DR/3000, DR 900 and DR/800 Series
Nume metoda:	SPADNS