



## Prepared agar plates m-ColiBlue24, pk/15

Nr. produs: 2805215

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

### m-ColiBlue24 Broth

m-ColiBlue24 Broth, Agar Plates

m-ColiBlue24 Broth, Agar Plates

Hach's ready-to-use media eliminates measuring, mixing, and autoclaving steps necessary to prepare dehydrated media. In 24 hours, Hach's m-ColiBlue24® Broth give you simultaneous results for total coliforms and E. coli. Glass ampules (26084-20), plastic ampules (260845-00), 100mL bottles (26084-42) of m-ColiBlue24® and agar plates (28052-15), afford maximum shelf life and ease of use for pre-measured and prepared media. Enumerate total coliforms and E. coli on one petri dish. Conduct emergency testing and routine monitoring. Get superior sensitivity - 1 CFU/100mL. Differentiate colonies easily - red and blue indicate total coliforms and blue specifies E. coli. Get optimal recovery of stressed and injured organisms. Minimize background growth of non-coliform bacteria. Eliminate special equipment and ultraviolet lamp. With ready-to-use m-ColiBlue24® Broth, you get these benefits:

#### 24-Hour Detection and Identification

Hach's m-ColiBlue24 Broth simultaneously detects and identifies both total coliforms and Escherichia coli (E. coli). Incubation time is 24 hours for complete enumeration. There is no confirmation step.

#### Easy to Read

m-ColiBlue24 makes it easy to differentiate between coliforms and E. coli—E. coli are blue; other coliforms are red; total coliforms are the sum of the two. It's that simple.

#### Easy to Interpret

m-ColiBlue24 requires no fluorescent lamps or color comparators—interpretation of the results is through easy, visual identification.

#### More Accurate

Enumerate test results directly with m-ColiBlue24 Broth. m-ColiBlue24 is much more accurate and precise than statistical-based tests such as Most Probable Number (MPN) and Multiple Tube Fermentation (MTF).<sup>1,2,3</sup>

#### Less Expensive

m-ColiBlue24 is the best value for an EPA-approved, 24-hour media for enumeration of total coliforms and E. coli. In addition, Hach Company not only offers the total solution, but also special discount packages on media, Petri dishes, membranes, and filtration equipment.

---

### Specificații

|             |            |
|-------------|------------|
| Cantitate:  | 15 /pk     |
| Material:   | Plastic    |
| Media Type: | m-ColiBlue |