



Plastic ampules, pseudomonas broth, pk/50

Nr. produs: 2812250

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

Pseudomonas broth plastic ampules, 2 mL. For detection of Pseudomonas species. pk/50.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Specificații

Cantitate:	50 50/pk
Material:	Plastic
Media Type:	Pseudomonas Broth