



Vacuum filter microfunnel, pk/50.

Nr. produs: 2831500

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Feature a completely assembled filtration unit: Funnel, filter, nutrient pad and petri dish. Filter: GN-6 Metrical with 0.45 μ m pore size and 47 mm diameter; white, gridded surface. Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Nutrient Pad: Cellulose. Fits 2831500, 2846100 and 2846150 manifolds. Sterile and individually bagged.

Specificații

Cantitate:	50
Diametru:	47 mm
Dimensiune por:	0.45 μ m
Material:	Mixed-Cellulose Ester