



Be Right™



Quality control standard inorganics for drinking water, 500 mL

Nr. produs: 2833049

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

Inorganics Quality Control Standard for Drinking Water. Ready to use. No dilution required. Mixed standard includes: Fluoride: 1 mg/L F^- , Nitrate: 2 mg/L NO_3^- -N, Phosphate: 2 mg/L PO_4^{3-} , Sulfate: 50 mg/L SO_4^{2-} . 500 mL Bottle.

Specificații

Concentrație:	Sulfate 50 mg/L SO_4^{2-}
Descriere:	Drinking Water Inorganics, Mixed Parameter QC Standard
Platforma:	Solu#ie