



Quality control standard metals for drinking water, low range, 500 mL

Nr. produs: 2833749

RON Preț (fără TVA): Contact

Disponibil

Ready to use - no dilutions required. Mixed standard includes: Copper: 1 mg/l Cu, Iron: 0.3 mg/L Fe, Manganese: 0.1 mg/L Mn. 500 mL Bottle.

Specificații

Concentrație:	Mangan: 0,1 mg/L Mn
Descriere:	Metals Drinking Water (Low Range) Mixed Parameter QC Standard
Platforma:	Solu#ie