



Sodium Standard Solution, 10 mg/L Na, 1 L

Nr. produs: 2835153

RON Preț (fără TVA): Contact

Disponibil

NaCl Standard, 10 mg/L Na, 1 L. For Polymetron 9240 & 9245 analyser.

Specificații

Cantitate:	1 L 1 L
Concentrație:	10 mg/L
Condiții de depozitare:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Instrument:	9245 & 9240
Parametru:	Sodium
Platforma:	Solu#ie