



Test kit color disc chlorine dioxide, model CLD-2

Nr. produs:

2842800

RON Preț (fără TVA):

Contact

Expediere în 2 săptămâni

Produse periculoase



Monitor CIO₂ levels on-site with the NEW Chlorine Dioxide Test Kit. The Chlorine Dioxide Test Kit uses the reliable and EPA-accepted DPD method to monitor CIO₂ concentration, while eliminating chlorine interference through the addition of glycine. With 100 tests per kit and available range of 0-2 mg/L and resolution of 0.1 mg/L, operators can accurately monitor around the 0.8 mg/L regulated limit.

Simple, cost-effective method for maintaining optimum chlorine dioxide levels in your plant

USEPA-accepted for waste water and potable water

Approved by some states for potable water under specific conditions (check with your regulatory agency)

Specificații

Domeniu de măsurare:	0 - 2.0 ppm ClO ₂
Interval de măsurare standard:	0 - 2.0 ppm ClO ₂
Model:	CLD-2
Model test kit:	CLD-2 ²
Numar teste:	100
Nume metoda:	Color Disc/DPD
Parameter:	Chlorine Dioxide
Parametru:	Chlorine Dioxide
Platforma:	Disc culoare
Smallest increment:	0.1
Tipul carcasei:	D