



Filter membrane, EZ-Pak, 0.45 micron, 47 mm, white, gridded, 600/pk (millipore)

Nr. produs: 2916800

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Millipore EZ Pack Membranes, 0.45 μ m pore size, 47mm

Millipore EZ Pack Membranes, 0.45 μ m pore size, 47mm

Dispenser provides superior convenience and assures sterility. Unique one press on the lever and a sterile membrane is presented (horizontally and inverted) for efficient and convenient assembly. The EZ Pak membranes come as a band of membranes that easily loads into the dispenser. Packaging materials are wound into the dispenser ensuring a clear operating area and quick change of bands. The membranes are Type HA, which signifies compliance with Standard Methods criteria.

Specificații

Cantitate:	600
Diametru:	47 mm mm
Dimensiune por:	0.45 μ m
Material:	Mixed-Cellulose Ester