



Bottle, wash, polyethylene, wide mouth, 1000 mL, 3/pk

Nr. produs: 2927204

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Polyethylene bottle with polypropylene closure and stem molded in one piece for leak-proof service. Tip can be cutback to increase flow. Available in Narrow or Wide Mouth.

Specificații

Cantitate:	3 pk
Capacitate:	1000 mL
Culoare:	Incolor
Material:	LDPE
Stem:	Transparent Cap/Stem
Unitate:	3 pk