



Bottle, Storage, Wide Mouth, Write-On, 1000 mL, 6/pk

Nr. produs: 2930004

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Large dedicated labeling areas help to standardize the way bottle contents are identified. A separate area is provided for name and a larger area allows for detail of contents or safety information and includes a defined space for date information. Bottles are molded of rugged HDPE in 500 and 1000 mL sizes and are available as narrow or wide mouth with leak proof polypropylene cap.

Specificații

Body Material:	High Density Plastic
Cantitate:	5 /bag
Capacitate:	1000 mL
Material:	HDPE
Package size UOM:	6pk