



Bottle BOD, 300 mL, disposable

Nr. produs:	2943101
RON Preț (fără TVA):	Contact
Expediere în 2 săptămâni	

Single 300 mL disposable BOD bottle without stopper. Plastic bottle is EPA approved for BOD measurement.

Box of 100 bottles available under article number 2943100. Plastic stoppers (25 pack) available under article number 2943900. Glass stoppers (quantity of 1) available under article number 190901.

Specificații

Body Material:	Plastic
Cantitate:	1 each
Capacitate:	300 mL
Material:	PET