



Sodium/Potassium Ionic Strength Adjustor (ISA) Powder Pillows, pk/100

Nr. produs:

4451569

Produse periculoase

RON Preț (fără TVA):
Expediere în 2 săptămâni

Contact



Ionic Strength Adjustor (ISA) for Sodium/Potassium determinations by Ion Selective Electrode (ISE) method. Powder form prevent from diluting samples. Ideal for low Sodium/Potassium level readings. Eliminate most interferences. Pack of 100 powder pillows.

Specificații

Cantitate:	100 /pk
Concentratie:	ISA
Package settings:	Powder Pillows
Parametru:	Sodium