



Be Right™



Beaker, griffin, low form, glass, 400 mL

| | |
|--------------------------|---------|
| Nr. produs: | 50048 |
| RON Preț (fără TVA): | Contact |
| Expediere în 2 săptămâni | |

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Specificații

| | |
|----------------|--------------------------|
| Capacitate: | 400 mL |
| Material: | Glass, Griffin, Low Form |
| Material corp: | Glass |
| Unitate: | 1 each |