



QL filters, pore size 8 - 12 μm , diameter 55 mm, circular, 100/ pk, qualitative

Nr. produs: 50652

ARTICOL ÎNVECHIT
Acest articol nu mai este disponibil.

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Specificații

Cantitate:	100 /pk
Diametru:	55 mm
Material:	Cellulose