



QL filters, pore size 8 - 12 μ m, diameter 125 mm, circular, 100/pk, qualitative

Nr. produs: 50657
RON Preț (fără TVA): Contact
Nu este afișată data expirării

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Specificații

Cantitate:	100 /pk
Diametru:	125 mm
Material:	Cellulose