



Cylinder, graduated, 5 mL $\pm$ 0.1 mL, 0.1 mL divisions (white markings)

Nr. produs: 50837
RON Preț (fără TVA): Contact
Disponibil

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes $\geq$ 25 mL.

Specificații

Cantitate:	1 each
Capacitate:	5 mL
Certificări:	Class B
Înălțime:	125 mm
Material corp:	Glass
Metoda de măsurare:	To Deliver
Tip:	Simplu Metric