



Cylinder, graduated, 25 mL +- 0.3 mL, 0.5 mL divisions (white markings)

Nr. produs: 50840
RON Preț (fără TVA): Contact
Expediere în 2 săptămâni

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Specificații

Cantitate:	1 each
Capacitate:	25 mL
Certificări:	Class B
Material corp:	Glass
Metoda de măsurare:	To Contain
Tip:	Simple Metric