



## Portable HQD rugged field case for rugged probes w/ extended cable length

Nr. produs: 8505600

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

### Rugged, Convenient Storage for Field Use

Rugged Field Case features custom cut foam inserts to hold meter, probes, and accessories in place. Holds portable meter and rugged probe with up to 30 meters of cable.

**All-in-One durable rugged field cases provide convenience and include space for sample collection containers, meter and rugged probes**

**Durable, die-cut foam inserts hold meter and probes securely in place**

---

### Specificații

Dimensiuni (Înălțime x Lățime x Adâncime): 420 mm x 458 mm x 216 mm