



EZ3508 Fluoride Analyser

Nr. produs: **EZ3508.XXXXXXXXXX**
RON Preț (fără TVA): Contact

Online ion-selective Fluoride measurements with standard addition for industrial applications

ISE technology for complex water matrices

The EZ3500 Series are used for water monitoring applications where ion-selective electrodes are the preferred analytical technique, but the water matrix is complex. Standard addition helps minimise matrix interferences that otherwise might contribute to the output signal of the electrode.

Discontinuous, automatic standard addition

The EZ3500 Series run discontinuous measurements with careful addition of a known standard in an analysis vessel. Discontinuous ISE analysis enhances control over the potentiometry, eliminates cross-contamination between cycles and reduces overall reagent consumption.

The EZ3500 Series combine unique technology with a set of analysis, control and communication features in an industrial analyser mainframe designed for the highest performance:

- Smart automatic features
- Control and communication via industrial panel PC
- Analogue and digital output options
- Multiple stream analysis (up to 8 streams)

There are many additional options available. Please contact Hach for more details.

Specificații

| | |
|-----------------------|---|
| Alarma: | 1x malfunctioning, 4x user-configurable, max. 24 VDC/0.5 A, potential free contacts |
| Alimentare electrică: | 100 - 240 VAC, 50/60 Hz |
| Calibrare: | Automatic, 2-point; frequency freely programmable |
| Calitate proba: | Maximum particle size 100 µm, < 0.1 g/L; Turbidity < 50 NTU |
| Certificări: | CE compliant / UL certified |
| Clasa de protecție: | Protection class: Analyser cabinet: IP44 / Panel PC: IP65 |
| Conținutul cutiei: | EZ3508 Fluoride Analyser, Instruction Manual, 1 x Double Bit Door Key, 1 x Mounting Brackets, 1 x ISE Electrode, 1 x empty 5L Reagent Container with Fittings (for F Standard Solution) and 1 x empty 10L Reagent Container with Fittings (for Buffer Solution) |

| | |
|--|--|
| Curatare automata: | Da |
| Cuve: | 8 minutes |
| Debitul probei: | 100 - 300 mL/min |
| Demineralised water: | For rinsing |
| Digital outputs: | Optional: Modbus (TCP/IP, RS485) |
| Dimensiuni (Înălțime x Lățime x Adâncime): | 690 mm x 465 mm x 330 mm |
| Domeniu de măsurare: | 5 - 100 mg/L F ⁻ |
| | Optional: |
| | 1 - 25 mg/L F ⁻ |
| | 2.5 - 50 mg/L F ⁻ |
| Earth connection: | Dry and clean earth pole with low impedance (< 1 Ohm) using an earth cable of > 2.5 mm ² |
| Esantion: presiune: | By external overflow vessel |
| Garanție: | 24 luni |
| Greutate: | 25 kg |
| Iesiri analogice: | Active 4 - 20 mA max. 500 Ohm load, standard 1, max. 8 (option) |
| Instrument air: | Dry and oil free according to ISA-S7.0.01-1996 quality standard for instrument air |
| Interferente: | Metal ions like aluminium > 72 mg/L, calcium > 108 mg/L and iron > 150 mg/L. Fats, oil, proteins, surfactants and tar. |
| Limita de detecție: | ≤ 1 mg/L |
| Material: | Hinged part: Thermoform ABS, door: PMMA |
| | Wall section: Galvanised steel, powder coated |
| Metoda de măsurare: | Discontinuous measurement by combined ion-selective electrode with standard addition |
| Model: | EZ3508 |
| Nr. flux probe : | 1, 2, 4 or 8 |
| Parametru: | Fluoride |
| Precizie: | Better than 2% full scale range for standard test solutions |
| Putere: | 100 - 240 VAC, 50/60 Hz |
| | Max. power consumption: 120 VA |
| Reagent Requirements: | Keep between 10 - 30 °C |
| Scurgere: | Atmospheric pressure, vented, min. Ø 32 mm |
| Temperatura ambiantă: | 10 - 30 °C ±4 °C deviation at 5 - 95% relative humidity (non-condensing) |
| Temperatura probă: | 10 - 30 °C |
| Validation: | Automatic; frequency freely programmable |

Conținutul cutiei

EZ3508 Fluoride Analyser, Instruction Manual, 1 x Double Bit Door Key, 1 x Mounting Brackets, 1 x ISE Electrode, 1 x empty 5L Reagent Container with Fittings (for F Standard Solution) and 1 x empty 10L Reagent Container with Fittings (for Buffer Solution)