



## Barium Chloride solution, 0.1 M, for EZ1036 Sulfate Analyser

Nr. produs:

APPC1036-02-EU

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Specificații

Cantitate: 2 L

Parametru: Sulfate