



Cylinder, Graduated, Mixing, Glass, 50 mL $\pm$ 0.4 mL, 1.0 mL divisions, glass stopper #16

Nr. produs: 104592

RON Preț (fără TVA): Contact

Low form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Low form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

Specificații

Cantitate:	1 each
Capacitate:	50 mL
Certificări:	Class B
Material corp:	Glass
Metoda de măsurare:	To Contain
Tip:	Simple Metric, Mixing