



Be Right™



General Boiler Water Test Kit, Model AL-94 MG-L (Drop Count Titration)

Nr. produs:

222403-EU

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Used for controlling boiler and cooling water treatment and checking corrosion, this test kit is an economical solution for the small boiler operator who needs to monitor alkalinity, chloride, hardness, iron, and pH.

Uses simple, economical drop count titration and color disc methods

Kit contains all required reagents and apparatus to perform analysis in the field

Uses drop count titration and color disc methods

Specificații

Aplicatie:	Boiler Water
Model:	AL-94 MG-L
Parametru:	Alkalinity, chloride, hardness, iron, and pH (4-10)
Platforma:	Disc culoare