



Hardness (Calcium and Total) Reagent Set, Digital Titrator

Nr. produs:

2272100-EU

Produse periculoase

RON Preț (fără TVA):
Disponibil

Contact



For Calcium and Total Hardness determination by EDTA Digital Titrator titration. Ranges: 10-160 and 100-4000 mg/L CaCO₃. Reagent Set contains: Buffer Solution (Hardness 1), CalVer 2 Calcium Indicator Powder Pillows, ManVer 2 Hardness Indicator Powder Pillows, Potassium Hydroxide Solution (8 N), EDTA Titration Cartridges (1 each) (0.800 and 0.0800 M), and 5.25N sulfuric acid.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specificații

Caracteristici speciale:	Reagent set, Digital Titrator
Domeniu de măsurare:	10 - 160 mg/L & 100-4000 mg/L CaCO ₃
Nume metoda:	EDTA Titration
Parametru:	Hardness (Calcium and Total)

Accesorii necesare

- Mănuși de protecție mărimea L, albastre, nitril, fără pudră, 50 de bucăți (Item SM743L)