



Test kit, nitrogen, nitrate color disc, model NI-14

Nr. produs:

1416100-EU

RON Preț (fără TVA):

Contact

Nu este afișată data expedierii

Produse periculoase



Utilizing the simple color disc method, this kit measures Nitrogen, Nitrate as NO_3^- and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Cei mai mici pași de incrementare:	0,05
Domeniu de măsurare:	0 - 1 mg/L NO_3^- -N
Domeniu de măsurare :	0 - 10 mg/L NO_3^- -N
Interval de măsurare standard:	0 - 1 mg/L NO_3^- -N; 0 - 10 mg/L NO_3^- -N
Model:	NI-14
Model test kit:	NI-14**
Numar teste:	50
Nume metoda:	Color disc/ Cadmium reduction
Parameter:	Nitrogen, Nitrate
Parametru:	Nitrogen, Nitrate - As NO_3^- -N
Platforma:	Disc culoare
Smallest increment:	0,05 0,5
Tipul carcasei:	B