



Be Right™



Test Kit, chlorine (free and total), model CN-80

Nr. produs:

2129000-EU

Produse periculoase

RON Preț (fără TVA):

Contact

Nu este afișată data expirării



Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specificații

Domeniu de măsurare:	0 - 0.68 mg/L Cl ₂
Domeniu de masurare :	0 - 10 mg/L Cl ₂
Interval de măsurare standard:	0 - 0.68 mg/L Cl ₂ ; 0 - 3.4 mg/L Cl ₂ ; 0 - 10 mg/L Cl ₂
Model:	CN-80
Model test kit:	CN-80 ²
Numar teste:	100
Nume metoda:	Color Disc/DPD
Parameter:	Chlorine (Free and Total)
Parametru:	Chlorine, low range - As free & total Cl ₂
Platforma:	Disc culoare
Smallest increment:	0.02
	0.1
	0.5
Tipul carcasei:	D