



Be Right™



Nitrate test kit, model NI-11

Nr. produs: 146803-EU

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Produse periculoase



One of Hach's most popular test kits, it uses the simple, low-cost color disc method.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specificații

Domeniu de măsurare:	0 - 40 mg/L NO ₃ ⁻ -N
Interval de măsurare standard:	0 - 40 mg/L NO ₃ ⁻ -N
Model:	NI-11
Model test kit:	NI-11**
Numar teste:	100
Nume metoda:	Color disc/ Cadmium reduction
Parameter:	Nitrogen, Nitrate
Parametru:	Nitrogen, Nitrate - As NO ₃ ⁻ -N
Platforma:	Disc culoare
Smallest increment:	1.0
Tipul carcasei:	B